**USER MANUAL
MANUAL DE USUARIO
BENUTZERHANDBUCH
MANUEL D'UTILISATION
MANUALE UTENTE
ユーザーマニュアル**

**KMINA - ALL-TERRAIN SCOOTER
KMINA - SCOOTER TODOTERRENO
KMINA - GELÄNDEROLLER
KMINA - SCOOTER TOUT-TERRAIN
KMINA - SCOOTERTUTTOTERRENO
KMINA - 全地形対応スクーター**

K10073

14

	ENGLISH	ESPAÑOL	DEUTSCH	FRANÇAIS	ITALIANO	日本語
1.	Handlebar	Manillar	Lenker	Guidon	Manubrio	ハンドルバー
2.	Rear-view mirror	Espejo retrovisor	Rückspiegel	Rétrovisueur	Specchio retrovisore	バックミラー
3.	Control panel	Panel de control	Bedienfeld	Panneau de commande	Pannello di controllo	コントロールパネル
4.	Cup holder	Portavasos	Getränkhalter	Porte-gobelets	Portabicchieri	カップホルダー
5.	Front lights	Luces delanteras	Vordere Lichter	Feux avant	Luci anteriori	フロントライト
6.	Front wheels	Ruedas delanteras	Vordere Räder	Roues avant	Ruote anteriori	フロントホイール
7.	Footrests	Reposapiés	Fußrasten	Repose-pieds	Poggiapiedi	フットレスト
8.	Rear wheels	Ruedas traseras	Hintere Räder	Roues arrière	Ruote posteriori	後輪
9.	Battery	Batería	Batterie	Batterie	Batteria	バッテリー
10.	Adjustable seat	Asiento ajustable	Verstellbarer Sitz	Siège réglable	Sedile regolabile	調整可能シート
11.	Charging port	Puerto de carga	Anschluss zum Aufladen	Port de charge	Porta di ricarica	充電ポート
12.	Folding armrest	Reposabrazos abatible	Klappbare Armlehne	Accoudoir rabattable	Bracciolo pieghevole	折りたたみ式アームレスト
13.	Folding backrest	Respaldo abatible	Umklappbare Rückenlehne	Dossier rabattable	Schienale pieghevole	折りたたみ式背もたれ
14.	Charger	Cargador	Ladegerät	Chargeur	Caricabatterie	充電器

WIDTH	HEIGHT	LENGTH
64 cm (25")	130 cm (51")	128 cm (47.5")
		

SEAT HEIGHT	SEAT DIMENSION	BACKREST HEIGHT	Ø REAR WHEEL	Ø FRONT WHEEL
49 cm (19")	50 x 48 cm (19.5"x19")	66 cm (26")	30 cm (12")	30 cm (12")

MAX. USER WEIGHT	WEIGHT WITHOUT BATTERY	WEIGHT WITH BATTERY
160 kg (350 lb)	79.5 kg (175 lb)	90 kg (198 lb)

EN | The KMINA - All-terrain Scooter has been designed to offer comfort and safety to the user when moving around. It supports a maximum weight of 350 lb and is designed for independent patients, elderly people or people with reduced mobility.

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EN | Press the lever on the handlebar structure to lift it and start assembly.

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B - SEAT / ASIENTO / SITZ / SIÈGE / SEDE / シート

EN | Insert the seat fully straight into the outer tube of the scooter until they are joined.

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EN | Press the lever on the handlebar structure to move it forward or backward and adjust according to user preference.

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


* **EN** | When the READY symbol flashes, this indicates a malfunction of the scooter (see annex for possible faults and solutions).

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EN | The brake lever allows you to switch between electric and manual driving modes.

■ **Electric driving (D):** When the brake lever is in the rear, in the electric driving position, the scooter remains in a braking state and cannot be pushed. At this time, the scooter can be driven electrically in a normal way.

■ **Manual driving (N):** When the brake lever is forward, in the manual driving position, the scooter is in brake release state and can be pushed. This function can only be activated when the scooter is on and off. During this time, you cannot drive the scooter.

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EN | Before driving, familiarise yourself with local pedestrian traffic regulations and practice. Make sure the brakes are working properly, all functions are normal and the battery is fully charged before using it on the road.

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EN | 1. Rotate the seat to get on the scooter:

- Pull up on the seat lock lever under the seat to release the seat.
 - Rotate the seat in the desired direction.
 - Release the seat lever and the seat will lock automatically.
2. Get on the scooter, turn the power lock key and switch the scooter on.
3. Adjust the speed regulator according to road conditions, environmental conditions and personal preference.
4. Hold the handlebar with both hands and slowly move the driving lever forward with your right fingers to advance the scooter.
5. When reversing, pay attention to the surrounding conditions. Use the fingers of your left hand to slowly move the drive lever forward. When parking, release the lever to bring the scooter to a smooth stop.
6. To stop the scooter, release the drive lever to return the scooter to the centre position and stop smoothly.

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Steel and plastic (ABS).

WARNINGS

- Make sure you are aware of the current legislation for the use of the scooter in your area or zone.
- DO NOT expose the scooter's electrical system to damp environments for long periods of time (rain, snow).
- DO NOT use in adverse weather conditions such as rain, snow, fog, high winds or extreme temperatures.
- DO NOT drive under the influence of alcohol, drugs, excessive fatigue or mental confusion.
- Avoid driving on terrain such as beaches, muddy areas or slippery surfaces.
- The scooter is designed for one person only. Do not carry passengers or use it to transport or tow objects.
- The rider should remain seated with both hands on the handlebars. Avoid sharp turns or sudden weight shifts to prevent tipping.
- Reduce speed in rough terrain or congested areas.
- Always turn off the ignition switch after parking and remove the key.
- When descending slopes, do not switch the brake lever to the manual drive position to enter manual driving mode, as there will be no braking effect. This is to ensure the safety of the user and the mobility scooter.
- If you detect faults, take the scooter to a professional technician.
- Read the safety instructions carefully and fully understand their contents before using the scooter.
- Reduce your speed in uneven terrain or congested areas.
- Always turn off the ignition switch after parking and remove the key.
- If you do not use the scooter for more than a week, disconnect the battery.
- When descending slopes, do not change the brake lever while using the scooter.
- DO NOT attempt to open, repair or modify the internal structure of the motor.
- DO NOT use gasoline or other chemicals to clean the housing or other parts of the electric scooter.
- Check daily to see if the battery needs charging.

- Inspect the brakes before riding to ensure that they are working properly.
- Regularly check the lever fixing bolts for looseness and ensure that the lever is in the neutral position.
- Avoid getting the electrical system wet, use the scooter only after it is completely dry.
- Periodically check the condition of the battery and clean the terminals if necessary.
- If the scooter is not used for an extended period of time, charge the battery approximately once a month.
- Regularly inspect the insulation of the cables and ensure that all connectors make good contact.
- Check the front and rear tyres for oscillation.
- Regularly inspect the condition of the scooter, such as brakes, battery, tyres and electrical connections.
- DO NOT cross puddles deeper than 2.5".
- Your surroundings may disturb the operation of devices that emit electromagnetic fields (shop alarm systems, automatic doors, etc.).
- The braking distance on slopes may be significantly greater than at ground level.
- If lit cigarettes are thrown onto a padded seating system, a fire may result causing damage, serious injury or death.
- Occupants of the electric vehicle are at particular risk of death or serious injury from such fires and the gases they produce, as they may not be able to move away from the electric vehicle.
- DO NOT smoke while operating this electric vehicle.
- Improper supervision or maintenance could result in injury, damage, or death due to ingestion or suffocation caused by parts or materials.
- Make sure you have all your belongings checked before turning on and using the scooter.
- Improper use of this product may result in injury or damage.
- DO NOT attach or tie a leash to the scooter to walk your pet.
- If you have any questions regarding warnings, cautions, or instructions, contact a healthcare professional or your dealer before attempting to use this equipment.
- Turn off the electric vehicle before getting on, off, or handling unwieldy objects.
- NEVER leave the mobility device on a slope.
- If you must brake in an emergency, simply release the drive lever and allow the mobility device to come to a complete stop.
- DO NOT switch off the electric vehicle in motion unless there is no other option.
- NEVER transport the electric vehicle with the occupant seated in the vehicle.
- Avoid storing or using the electric vehicle near fire or combustible products.
- DO NOT exceed the maximum permissible load.
- When servicing, repairing or lifting any part of the electric vehicle, be aware of the weight of the individual components, especially the batteries. Ensure that you always adopt the correct lifting position and ask for assistance if necessary.
- Ensure that moving parts of the mobility device, such as wheels or the seat lift, do not cause injury, especially when children or animals are nearby.
- The scooter is NOT recommended for use by minors.
- DO NOT expose the mobility scooter to direct sunlight for long periods of time. Metal surfaces and parts such as the seat may overheat.
- DO NOT connect any electrical device to the electric vehicle that has not been expressly certified by IMPROVING MOBILITY DEVELOPMENT S.L. for this purpose.
- Improper use could cause the electric vehicle to smoke, spark or catch fire. Damage, serious injury or death could result from fire.
- DO NOT use the electric vehicle for any purpose other than that for which it was intended.
- If the electric vehicle starts to smoke, spark or burn, stop using it and have it repaired immediately.
- To avoid electric shock, check the plug and cord for cuts or frayed wires. Replace cut or frayed cords immediately.
- Corroded electrical components due to exposure to water or other liquids can cause damage, serious injury or death.
- Reduce exposure of electrical components to water or other liquids.
- Electrical components damaged by corrosion should be replaced immediately.
- Electrical vehicles that are frequently exposed to water or other liquids may require more frequent replacement of electrical components.
- NEVER cover the lighting system with fabrics.
- DO NOT insert fingers between scooter parts. The use of gloves is recommended when handling scooter parts.
- Exercise extreme caution with metal and heavy parts.
- A fault in the electrical system may cause unusual behaviour, such as the light coming on steady or not coming on at all.
- Only climb gradients that do not exceed the maximum safety slope and only with the backrest in an upright position.
- Avoid braking sharply or accelerating on slopes.
- NEVER attempt to negotiate an obstacle when going up or down a slope.
- Always approach obstacles in a straight line. Make sure that the front and rear wheels go over the obstacle in one go, do not stop right in the middle.
- Avoid varying your centre of gravity as well as sudden changes of direction when the mobility device is in motion.

- Always distribute the weight evenly when carrying the mobility device. Always try to keep the centre of gravity of the electric vehicle in the middle and as close to the ground as possible.
- Be aware that the electric vehicle will brake or accelerate if the driving speed changes while in motion.
- Drive in narrow spaces at minimum speed and with due caution.
- There is an increased risk of rollover when making turns.
- Reduce speed before making turns. Only accelerate when you have completed the turn.
- Please note that the seat height has a major influence on the centre of gravity.
- The higher the seat height, the greater the risk of tipping over.
- Drive with extreme caution on soft surfaces, especially on uphill and downhill slopes.
- During the process, pay particular attention to the tipping stability of the electric vehicle.
- DO NOT attempt to carry out maintenance work that is not described in this user manual. Such repair and/or maintenance must be carried out by a qualified technician. Contact a dealer.
- For safety reasons and to prevent accidents resulting from inadvertent wear and tear, it is important that this electric vehicle is serviced annually under normal operating conditions.
- Under difficult operating conditions, e.g. daily driving on steep slopes or in the case of use in medical care with frequent changes of wheelchair users, it is advisable to carry out intermediate checks on the brakes, accessories and running gear.
- Spare parts must match the original parts.
- Always state the serial number of the electric vehicle to make it easier to order the correct spare parts.
- Only use batteries that have been approved by IMPROVING MOBILITY DEVELOPMENT S.L. for this electric vehicle.
- Continued use of an electric vehicle that has not been adjusted to the correct specifications could cause erratic vehicle behaviour resulting in damage, serious injury or death.
- After any adjustment, repair or maintenance and before proceeding with use, make sure that all fastening parts are assembled and tightened securely.
- Ensure that the handgrip is correctly adjusted before riding the scooter.
- After making adjustments to the angle of the handgrip and before use, the handgrip must be firmly locked in position. Otherwise, the scooter could fall, resulting in personal injury and/or damage to the vehicle. Gently push or pull on the handle to ensure that the handle is securely attached to the adjustment plate.
- NEVER drive around obstacles at a steep angle.
- Before overcoming obstacles, place the backrest in an upright position.
- DO NOT attempt to get off the scooter on a slope or incline.
- Protect the battery charger from moisture. Always charge the batteries in a dry environment.
- DO NOT use the battery charger if it has been dropped or damaged.
- NEVER attempt to recharge the batteries by connecting the cables directly to the battery terminals.
- Only use an extension cable if absolutely necessary. If you must use an extension cable, make sure it is in perfect condition.
- Under no circumstances may this electric vehicle be used as a vehicle seat or to transport the user in a vehicle.
- Risk of twisting when lifting heavy parts. Use the correct lifting techniques.
- When cleaning electronic components, observe their protection class with regard to water ingress.
- Make sure that no water splashes onto the plug or wall socket.
- DO NOT touch the socket with wet hands.
- All cleaning agents and disinfectants used must be effective, compatible with each other and must protect the materials to be cleaned.
- Check that all wheels, bearings, handlebars, seat and other parts are in good condition.
- NEVER use corrosive fluids (alkalis, acids, etc.) or abrasive cleaning agents.
- We recommend using a normal household cleaner, such as liquid soap, unless otherwise specified in the cleaning instructions.
- DO NOT use solvents (cellulose strippers, acetone, etc.) that change the structure of the plastic or dissolve the attached labels.
- Always make sure that the product has dried completely before using it again.
- When the scooter is not in use, it is recommended that it be kept indoors or protected outdoors with a cover and properly supervised.
- The product must not be cleaned in automatic washing systems, using steam or pressure cleaning equipment.
- Dirt, sand and salt water may damage the bearings and steel parts may rust if the surface is damaged.
- This product may contain substances that could be harmful to the environment if it is disposed of in places (landfills) that are not legally permitted.
- DO NOT dispose of the battery with household waste.
- DO NOT dispose of the battery in a fire.
- The battery must be taken to a suitable collection point. There is no charge for returning the battery and it is required by law.
- Only dispose of the discharged battery.
- Cover the lithium battery with non-flammable packing material before disposal.

- The scooter has an interference resistance rating of 20 v/m, which would provide useful protection against common sources of radiated EMI.
- DO NOT use a portable CB transceiver or communication device such as a cellular phone when the scooter is powered on.
- Be careful and keep away from electromagnetic wave releasing equipment such as radio, TV station, etc.
- If unintentional movement occurs or the brake is released, turn off the motorized scooter as soon as it is safe.
- Note that adding accessories or components, or repositioning the scooter may make it more susceptible to interference from radio wave sources.
- Report all incidents of unintended movement or brake release to the scooter manufacturer and note whether there is a radio wave source nearby.
- Check the condition of the battery to decide whether it needs to be charged.
- Check the wheel fixing bolts and nuts regularly for looseness.
- Check tyre pressure and degree of wear regularly to see if tyres need to be replaced.
- Check the seat for looseness or sudden shaking before driving.
- Check whether the brake meets the braking requirements and the scooter can be placed on the appropriate slope for further inspection.
- Check the motor to see if it needs to be replaced, especially after using it for one year.
- Clean the surface of the scooter frequently and do not wash the scooter directly with water to avoid permanent damage. Wipe the scooter with a damp cloth.
- Develop good driving habits and DO NOT operate the scooter beyond the permitted range.
- The scooter should be stored in a dry and ventilated place. Please do not place it in water or in a very humid environment (humidity should be less than 85%).
- If the scooter has not been used for a long time, take dust prevention measures.
- The air tyre pressure should be checked at least once a week.
- Proper pressures will prolong the life of your tyres and ensure smooth operation when driving the scooter. It is vitally important to maintain tyre pressure at 30-25 psi (2-2.4 bar) at all times.
- The KMINA scooter does not include a tyre inflator.

Some parts of the scooter are susceptible to temperature changes. The control panel can only operate safely between -25 °C and 50 °C.

- The batteries may freeze at extremely low temperatures and the scooter may not be able to operate at such low temperatures. If it is at extremely high temperatures, it may operate at a lower speed due to a safety feature of the controller that prevents damage to the motor and other electrical components.
- IMPROVING MOBILITY DEVELOPMENT S.L. is not liable for its own damage, damage to third parties or damage to property, whether accidental or in bad faith.
- IMPROVING MOBILITY DEVELOPMENT S.L. is not liable for theft or loss of parts or the scooter.
- The CE marking will be invalidated if components or accessories are replaced or added that have not been approved by IMPROVING MOBILITY DEVELOPMENT S.L. for this product.
- Under no circumstances should you modify, add, remove or damage any part or function of the mobility scooter, otherwise personal injury and damage to the mobility scooter may result.
- Use the original product accessories supplied by us, untested or inapplicable accessories may cause failure.
- The company adding or replacing components or accessories to the product is responsible for the conformity assessment/ CE marking or for the registration of the electric scooter as a special design, as well as for the related documentation.
- NOT suitable for transport by plane.

CHARGING AND MAINTENANCE OF THE BATTERY

- Read all instructions carefully before using the battery charger.
- DO NOT charge the scooter outdoors, in confined spaces or under high temperature, flammable gas, fog or steam.
- DO NOT place the charger on the seat during charging.
- Ensure the scooter is switched off before charging.
- Ensure that the charger plug is dry and undamaged before connecting it to the scooter.
- After charging, remember to turn the charging lid backwards to prevent it from getting wet.
- The charger should be kept dry in an environment between -13°F and 104°F.
- Only use our original charger. A charger of another brand may shorten the warranty period, and an unknown charger may cause damage and injury.
- In all cases, the charger must only be repaired by an authorised dealer.

MAINTENANCE AND CLEANING

- Check the tyre pressure and degree of wear regularly to see if the tyres need to be replaced.

- Check the seat if it is loose or shakes sharply before driving.
- Check whether the brake meets the braking requirements and the scooter can be placed on the appropriate slope for further inspection.
- Check the motor to see if it needs to be replaced, especially after using it for one year.
- Clean the surface of the scooter frequently and do not wash the scooter directly with water to avoid permanent damage. Wipe the scooter with a damp cloth.
- Develop good driving habits and operate the scooter beyond the permitted range.
- The scooter should be stored in a dry and ventilated place. Please do not place it in water or in a very humid environment (humidity should be less than 85%).
- If the scooter has not been used for a long time, take dust prevention measures. Keep the battery fully charged and charge it at least once a month.

PERFORMANCE PARAMETERS

- Maximum speed: 9,3 mi/h.
- Maximum slope: 12°.
- Braking method: Electromagnetic.
- Lithium battery: 12V/52Ah.
- Motor: 24V/800W.
- Charger: 24V/5A.
- Charger cable: 72".
- Control panel: 24V/120A.

WARRANTY

IMPROVING MOBILITY DEVELOPMENT S.L. guarantees that this device is free from defects.

- The warranty on the engine, control panel, battery and charger shall be valid for 1 year from the date of purchase of the product.
- The product warranty shall remain in effect for 3 years from the date of original consumer purchase. The warranty does not extend to parts exposed to natural wear and tear that require periodic replacement. This warranty does not cover:
 - Replacement parts supplied by anyone other than IMPROVING MOBILITY DEVELOPMENT S.L. or the authorised distributor.
 - Any failure of the unit during the warranty period that is not caused by a defect in material or workmanship or if the failure is caused by misuse.
 - IMPROVING MOBILITY DEVELOPMENT S.L. is not liable for its own damage, damage to third parties or damage to property, whether accidental or in bad faith.
 - IMPROVING MOBILITY DEVELOPMENT S.L. is not liable for theft or loss of parts or the scooter.

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EN | POSSIBLE FAILURES AND SOLUTIONS

FAILURES	POSSIBLE CAUSES	SOLUTIONS
The driving distance has drastically reduced	Batterie not charged for enough time	Charge batteries for eight hours or more
	Batterie is aging to low capacity	Replace batterie
Battery cannot charge or battery gauge shows empty after charging	Battery fault	Replace batterie
	Charge fault	Replace batterie
	Charger plug damaged	Check plug
	Loose connection	Try a wall socket in a different room
	No output from wall outlet	Unplug from wall & change fuse
	Fuse in charger mains plug blown	Switch off and press button back in
	Button on battery pack has popped out	Switch off and press button back in
	Output fuse in charger blown	Unplug from wall and contact dealer
Battery charging current is high	Battery fault	Replace batterie
	Scooter switched on during charging	Switch off the scooter
The scooter cannot work	Brake-release lever disengaged	Engage brake-release lever
	Flat battery	Charge batterie pack
	Scooter is not switched on with key switch	Ensure the key switch is switched on
	The batterie cable is loose	Check batterie cable
	Charger plugged in	Unplug charger
	Button on batterie pack popped out	Reset circuit-breaker button
	Disconnected loom or plugs	Check all plugs & looms
	Control system fault	Contact with dealer
Motor runs irregularly and/or noisily	Electrical malfunction	Contact with dealer
	Control system fault	Contact with dealer
<p>* Do not attempt to disassemble any parts of the scooter control system, batterie or batterie charger. * The control system is critically important and users should not deal with it arbitrarily.</p>		

The scooter is engineered with a self-diagnostic controller which continuously monitors the operating conditions of your scooter. If it detects a problem, it will indicate error message by flashing the READY on the display and sounding audible will beep simultaneously. You should count the number of flashes and beeps separated by a short delay between each sequence, and check the list below to see what kind of problem has happened according to the number.

NUMBER OF FLASHES/ BEEPS	FAILURES	POSSIBLE CAUSES	SOLUTIONS
1	Low battery	The battery is discharged	The battery needs to be charged
2	Low voltage of battery	Lack of power	The battery needs to be charged
3	High voltage of battery	Too high voltage or overload	Decrease the speed
			Check the battery connection
4	Electrical current above limit	Electric current above the limit of the motor	Check the connections of the motor and related wiring
			Turn off, wait a few minutes and then
5	Braking issue	Manual driving mode is activated	Confirm that the lever is in the correct position
6	Driving levers activated	When the scooter is switched on, the position of the driving levers is not static	Make sure that the driving lever is in the neutral position
			It may be necessary to the driving lever may need to be recalibrated
7	Several failures	Possible faults in the throttle, driving levers or their wiring	Check the condition of the throttle, drive levers and their wiring
8	Motor voltage problem	Problem with the motor and other cables	Check the motor and other cables
9	Others	Internal problems in the control panel	Check all connection and wiring
10	Sliding problems	The pushing speed or sliding is too limited	Turn the scooter on and off

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