Español

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English

INSTRUCTIONS FOR USE AND PRESERVATION

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

REGULATION

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthosis.

SUITABLE FOR

- Weakness or inability to dorsiflex the ankle (dropfoot) as a result of neurological injuries, congenital alterations or muscular disorders.
- Flaccid paralysis associated with hemiparesis.

FITTING INSTRUCTIONS

THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.

In order to obtain a greater degree of therapeutic effectiveness in the different pathologies and to prolong the life and utility of the product, it is fundamental that the correct size adequate for each patient or user is chosen (enclosed in the packet is a diagram with all the sizes and the equivalent in centimetres). An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

To fit observe the following aspects:

Using an interface (i.e. a sock) is recommended.

A- ADJUSTING THE PERIMETER

A.1-The product can be adapted to a specific ankle perimeter by adjusting the part containing the elastic guide (pin).

A.2-Depending on the patient's ankle perimeter in the malleolus area, position this part between the first and last of the 4 holes on the supramalleolar part. See the following table for guidance:

Size	• (cm)	• (cm)	(cm)	(cm)
1	14	16	18	20
2	20	22	24	26
3	26	28	30	32

B- SECURING THE HOOKS

B.1-If using the anti-slip part with a hook, until the laces on the footwear and position it between the shoe tongue and laces so it is centred. Secure the lower micro-hook by holding the laces.

B.2-If using the "S" hooks, you must insert them in the eyelets with the longest part inside the shoe. The elastic traction strap must go through both hooks.

C- ADAPTING THE SUPRAMALLEOLAR PART

C.1-Open the supramalleolar part until you hear a click; the part will be flat.

C.2-Position the supramalleolar part in the back of the ankle so it coincides with the centre of the silicone pad against the Achilles tendon.

C.3-Press the ends of the part toward the leg until you hear a click which unlocks the part to be wrapped around the ankle.

C.4-Secure the tip of the widest micro-hook so the buckle is between the 2 narrowest straps.

C.5-Pull and secure the narrowest straps.

C.6-Insert the elastic traction band through the guide (pin) on the front.

C.7-Secure the elastic band to the hook (or hooks) and pull on the ends with the micro-hook and then secure them to get the desired traction.

D- FITTING THE INNER STRIP Ref.: AB12/AB22

D.1-Secure the strip in the middle of the foot with the Velcro closures in the instep. Make sure the closure hook is in the right direction. Then, insert the elastic traction band through the guide on the front of the supramalleolar part. Secure it to the hook and pull on the ends to get the desired traction.

E- INSERTING THE SHIN SUPPORT Ref. AB14/AB24

E.1-The shin support must correspond to the size of the Boxia AB100, after being previously adapted and checked for the patient.

E.2-Secure using the micro-hook on the inside and front of the supramalleolar part.

E.3-The "Y"-shaped shin support has 3 ends. The shortest one along with the longest and narrowest will hug the shin at the height of the top of the fibula. The remaining end will go over the front of the tibia to be secured to the supramalleolar part on the inside.

PRECAUTIONS

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to sweat, we recommend the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat.

Products marked with the  symbol contain natural rubber latex and can cause allergic reactions in latex-sensitive individuals. Look at the back of the packaging for products marked with the symbol.

RECOMMENDATION-WARNINGS

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthosis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed.

For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions.

MANUFACTURE-CHARACTERISTICS

The material features used in production have been tested and approved and meet all the European quality requirements.

All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.

RECOMMENDATIONS FOR PRESERVATION AND CLEANING

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap. To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product. You can store it completely flat (fig. F.1) or in its circular form (fig. F.2).

**GUARANTEES**

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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